

Analyte	R _t (mins)	Adduct	Precursor Ion m/z	Quant Ion m/z	Conf Ion m/z	Analyte	R _t (mins)	Adduct	Precursor Ion m/z	Quant Ion m/z	Conf Ion m/z
Pirimicarb	4.24	[M+H] ⁺	239.2	182.1	72.0	Spiromesifen	8.66	[M+NH ₄] ⁺	388.1	273.1	187.0
Pirimicarb Desmethyl	2.71	[M+H] ⁺	225.1	168.2	72.1	Spiromesifen Alcohol	5.01	[M+H] ⁺	273.2	187.1	179.1
Pirimiphos Methyl	7.34	[M+H] ⁺	306.1	164.2	108.1	Spirotetramat	6.38	[M+H] ⁺	374.2	302.3	216.2
Prallethrin	7.69	[M+H] ⁺	301.2	133.0	151.2	Spiroxamine	5.95	[M+H] ⁺	298.3	144.2	100.2
Prochloraz	7.39	[M+H] ⁺	376.0	308.1	70.1	Sulfoxaflor	2.39	[M+NH ₄] ⁺	295.2	174.1	154.1
Profoxydim	7.71, 9.00	[M+H] ⁺	466.2	280.0	180.0	Sulprofos	8.56	[M+H] ⁺	323.0	219.1	139.1
Promecarb	5.88	[M+H] ⁺	208.1	109.0	151.1	TCMTB	5.48	[M+H] ⁺	239.0	180.0	136.0
Propamocarb	1.41	[M+H] ⁺	189.1	102.0	144.0	Tebufenozide	6.78	[M+H] ⁺	353.2	133.0	104.8
Propaquizafop	8.21	[M+H] ⁺	444.1	299.2	371.2	Tebufenpyrad	8.19	[M+H] ⁺	334.2	117.1	145.1
Propargite	8.74	[M+NH ₄] ⁺	368.2	231.2	175.1	Tebuthiuron	3.89	[M+H] ⁺	229.1	172.0	116.0
Propetamphos	6.13	[M+H] ⁺	282.1	138.1	156.1	Tepraloxymid	4.10, 6.19	[M+H] ⁺	342.2	250.1	166.1
Propoxur(S)	3.69	[M+H] ⁺	210.1	168.2	111.1	Terbufos Sulfone	5.46	[M+H] ⁺	321.0	115.0	143.0
Prosulfuron	5.29	[M+H] ⁺	420.1	167.1	141.1	Terbufos Sulfoxide	5.49	[M+H] ⁺	305.1	97.0	187.0
Pymetrozine	1.44	[M+H] ⁺	218.1	105.1	78.1	Terbutylazine	5.71	[M+H] ⁺	230.1	174.1	104.1
Pyraclostrobin	7.30	[M+H] ⁺	388.1	163.1	194.1	Tetrachlorvinphos	6.86	[M+2+H] ⁺	366.9	127.1	206.0
Pyraflufen Ethyl	7.13	[M+H] ⁺	413.0	339.0	253.1	Tetramethrin	7.91, 8.10	[M+H] ⁺	332.2	164.1	135.1
Pyrazophos	7.31	[M+H] ⁺	374.1	222.2	194.1	Thiabendazole	2.48	[M+H] ⁺	202.0	175.0	131.1
Pyridaben	9.22	[M+H] ⁺	365.1	309.0	147.1	Thiacloprid	2.55	[M+H] ⁺	253.0	126.1	99.1
Pyridalyl	10.21	[M+2+H] ⁺	492.0	110.9	164.0	Thiamethoxam	1.65	[M+H] ⁺	292.0	211.1	181.1
Pyrimethanil	5.45	[M+H] ⁺	200.1	107.1	168.1	Thifensulfuron Methyl	3.28	[M+H] ⁺	388.0	167.1	205.0
Pyriproxyfen	8.39	[M+H] ⁺	322.1	96.0	227.1	Thiobencarb	7.46	[M+H] ⁺	258.1	125.0	89.0
Quinalphos	6.78	[M+H] ⁺	299.1	163.1	147.1	Thiodicarb	4.34	[M+H] ⁺	355.1	163.2	88.1
Quinoxifen	8.50	[M+H] ⁺	308.0	197.1	214.1	Thionazin	4.74	[M+H] ⁺	249.1	193.1	97.0
Quizalofop Ethyl	8.01	[M+H] ⁺	373.1	299.2	255.1	Topramezone	1.63	[M+H] ⁺	364.1	334.1	125.1
Resmethrin	9.40	[M+H] ⁺	339.2	128.1	171.1	Triadimefon	6.07	[M+H] ⁺	294.1	197.0	225.0
Rimsulfuron	3.94	[M+H] ⁺	432.1	182.1	139.0	Triadimenol	6.25	[M+H] ⁺	296.1	70.2	99.0
Rotenone	6.71	[M+H] ⁺	395.2	213.2	192.1	Triazophos	6.19	[M+H] ⁺	314.1	162.1	119.1
Saflufenacil	5.32	[M+H] ⁺	501.1	349.1	198.0	Tribenuron Methyl	4.59	[M+H] ⁺	396.1	155.1	181.1
Sedaxane	6.20, 6.54	[M+H] ⁺	332.2	159.0	139.0	Trichlorfon	2.26	[M+H] ⁺	256.9	109.0	221.0
Sethoxydim	8.03	[M+H] ⁺	328.2	178.0	220.1	Tricyclazole	2.80	[M+H] ⁺	190.0	163.1	136.1
Simazine	3.66	[M+H] ⁺	202.1	104.1	132.1	Trifloxystrobin	7.78	[M+H] ⁺	409.1	186.2	206.2
Spinetoram	8.14	[M+H] ⁺	748.5	142.1	203.1	Triflumizole	7.87	[M+H] ⁺	346.1	278.0	73.0
Spinosad A	7.69	[M+H] ⁺	732.5	142.1	98.0	Triforine	5.23	[M+2+H] ⁺	434.9	213.0	98.2
Spinosad D	8.10	[M+H] ⁺	746.5	142.1	98.0	Zoxamide	7.09	[M+H] ⁺	336.0	187.0	159.0
Spirodiclofen	8.91	[M+H] ⁺	411.1	313.1	71.1						