

Analyte	R _t (mins)	Adduct	Precursor Ion m/z	Quant Ion m/z	Conf Ion m/z	Analyte	R _t (mins)	Adduct	Precursor Ion m/z	Quant Ion m/z	Conf Ion m/z
3-OH Carbofuran	2.25	[M+H] ⁺	238.1	181.2	163.1	Cyprosulfamide	3.30	[M+H] ⁺	375.1	135.1	254.1
5-OHThiabendazole	1.66	[M+H] ⁺	218.0	147.2	191.1	Cyromazine	1.15	[M+H] ⁺	167.1	125.2	68.2
Abamectin	9.45	[M+NH ₂] ⁺	890.5	305.3	567.5	DEF	9.20	[M+H] ⁺	315.1	169.0	113.0
Acephate	1.26	[M+H] ⁺	184.0	143.1	125.1	Demeton-S sulfone	2.55	[M+H] ⁺	291.1	235.1	263.1
Acetamiprid	2.24	[M+H] ⁺	223.1	126.1	90.1	Dialifos	7.46	[M+H] ⁺	394.0	208.1	181.0
Aldicarb	2.95	[M+NH ₂] ⁺	208.1	116.1	89.0	Diazinon	7.12	[M+H] ⁺	305.1	169.1	153.2
Aldicarb sulfone	1.44	[M+NH ₂] ⁺	240.1	148.0	86.0	Diazinon OA	5.32	[M+H] ⁺	289.1	153.2	233.1
Aldicarb sulfoxide	1.37	[M+NH ₂] ⁺	224.1	132.0	89.1	Dichlormid	3.85	[M+H] ⁺	208.0	140.0	81.2
Allethrin	8.33	[M+H] ⁺	303.2	135.1	123.1	Dichlorvos	3.63	[M+H] ⁺	221.0	109.1	127.0
Ametoctradin	7.64	[M+H] ⁺	276.2	149.1	176.2	Dicrotophos	1.87	[M+H] ⁺	238.1	112.2	193.1
Atrazine	4.64	[M+H] ⁺	216.1	174.0	104.0	Diethofencarb	5.53	[M+H] ⁺	268.2	124.1	180.2
Azinphos ethyl	6.30	[M+H] ⁺	346.0	132.1	223.0	Diflubenzuron	6.66	[M+H] ⁺	311.0	158.0	141.0
Azinphos methyl	5.14	[M+H] ⁺	318.0	132.0	124.9	Dimethenamid	5.70	[M+H] ⁺	276.1	244.1	168.2
Azinphos methyl OA	2.98	[M+H] ⁺	302.0	132.2	160.0	Dimethoate	2.23	[M+H] ⁺	230.1	199.0	125.0
Azoxystrobin	5.59	[M+H] ⁺	404.1	372.1	344.1	Dimethomorph	5.76, 6.07	[M+H] ⁺	388.1	301.0	165.1
Bendiocarb	3.72	[M+H] ⁺	224.1	167.1	109.1	Dinotefuran	1.36	[M+H] ⁺	203.1	129.1	114.2
Benoxacor	5.23	[M+H] ⁺	260.1	134.1	120.1	Dioxacarb	2.26	[M+H] ⁺	224.1	123.1	167.1
Bifenazate	6.27	[M+H] ⁺	301.1	198.0	170.1	Dioxathion	8.10	[M-C ₂ H ₄ O ₂ PS ₂] ⁺	271.1	97.0	125.0
Bitertanol	7.41	[M+H] ⁺	338.2	269.3	99.1	Disulfoton sulfone	4.59	[M+H] ⁺	307.0	261.1	125.0
Boscalid	5.85	[M+H] ⁺	343.0	307.0	140.0	Disulfoton sulfoxide	4.49	[M+H] ⁺	291.0	185.1	213.1
Bupirimate	6.68	[M+H] ⁺	317.2	210.2	237.3	Diuron	4.82	[M+H] ⁺	233.0	72.1	160.0
Buprofezin	8.24	[M+H] ⁺	306.1	201.1	106.1	DMST	3.90	[M+H] ⁺	215.1	106.1	151.0
Cadusafos	7.58	[M+H] ⁺	271.1	159.0	131.0	Dodine	7.56	[M+H] ⁺	228.3	186.3	60.1
Carbaryl	4.07	[M+NH ₂] ⁺	219.1	145.1	127.0	Emamectin	8.57	[M+H] ⁺	886.5	158.1	126.1
Carbendazim	2.10	[M+H] ⁺	192.1	160.1	132.1	Ethiofencarb	4.27	[M+H] ⁺	226.1	107.1	169.1
Carbofuran	3.77	[M+H] ⁺	222.1	165.2	123.2	Ethiofencarb sulfone	1.90	[M+NH ₂] ⁺	275.1	107.1	201.1
Carboxin	3.97	[M+H] ⁺	236.1	143.0	93.0	Ethiofencarb sulfoxide	1.98	[M+H] ⁺	242.1	107.1	185.0
Carfentrazone ethyl	6.88	[M+H] ⁺	412.0	346.1	366.0	Ethion	8.31	[M+H] ⁺	385.0	199.1	143.0
Chlorantraniliprole	5.24	[M+H] ⁺	484.0	286.0	194.0	Ethion monoxon	6.73	[M+H] ⁺	369.0	199.0	143.0
Chlorfenvinphos	7.21	[M+H] ⁺	359.0	170.0	99.1	Ethiprole	5.77	[M+NH ₂] ⁺	413.9	351.0	255.0
Chlorimuron ethyl	5.73	[M+H] ⁺	415.1	186.0	83.0	Ethofumesate	5.55	[M+H] ⁺	287.1	121.1	241.1
Chlorpyrifos	8.47	[M+H] ⁺	349.9	198.0	97.0	Ethoprop	6.46	[M+H] ⁺	243.1	173.0	131.0
Chlorpyrifos OA	6.65	[M+H] ⁺	334.0	278.0	197.9	Etofenprox	9.75	[M+NH ₂] ⁺	394.2	177.2	107.1
Clethodim	7.71	[M+H] ⁺	360.3	164.1	136.1	Etoxazole	8.73	[M+H] ⁺	360.2	141.0	304.2
Clofentezine	7.38	[M+H] ⁺	303.0	138.1	102.0	Famoxadone	7.24	[M+NH ₂] ⁺	392.2	331.1	238.0
Cloransulam methyl	4.13	[M+H] ⁺	430.0	398.1	370.0	Fenamidone	5.76	[M+H] ⁺	312.1	236.1	92.2
Clothianidin	1.99	[M+H] ⁺	250.0	169.1	132.0	Fenamiphos	6.71	[M+H] ⁺	304.1	217.1	202.0
Coumaphos	7.07	[M+H] ⁺	363.0	227.1	307.1	Fenamiphos sulfone	4.10	[M+H] ⁺	336.1	266.1	188.1
Crotoxyphos	5.86	[M+NH ₂] ⁺	332.1	127.1	193.1	Fenamiphos sulfoxide	3.96	[M+H] ⁺	320.1	233.1	171.1
Crufomate	6.77	[M+H] ⁺	292.1	236.1	108.1	Fenazaquin	9.21	[M+H] ⁺	307.2	161.2	57.2
Cyantraniliprole	4.33	[M+2+H] ⁺	475.0	286.0	444.1	Fenhexamid	6.39	[M+H] ⁺	302.1	178.0	97.2
Cyazofamid	6.52	[M+H] ⁺	325.1	108.1	261.2	Fenobucarb	5.49	[M+H] ⁺	208.1	95.0	152.0
Cyflufenamid	7.42	[M+H] ⁺	413.1	295.1	203.0	Fenoxaprop ethyl	8.04	[M+H] ⁺	362.1	288.1	91.1
Cymoxanil	2.48	[M+H] ⁺	199.1	128.1	111.1	Fenoxycarb	6.80	[M+H] ⁺	302.1	88.1	116.1
Cyphenothrin	9.27	[M+NH ₂] ⁺	393.2	151.2	123.2	Fenpropimorph	6.42	[M+H] ⁺	304.3	147.2	119.1

